



PRODUCT DATA SHEET

PRODUCT: Research Tool

CATALOG NUMBER: RT 01001-CLTH

QUANTITY & CONCENTRATION: Soluble in DMSO 25 mg/ml (20 µl; 100 mM); MW = 754.8

APPEARANCE: Yellowish-orange liquid

SHIPPED IN: Suspended in DMSO, ≤ -20°C

BACKGROUND/DESCRIPTION:

RT 01001-CLTH is a bispecific molecule that binds extracellular domain of EGFRvIII.

Importantly:

Maximum Excitation: 430 nm

Maximum Emission: 520 nm

Despite the fluorescent staining, due to its affinity towards specific receptors, prolonged incubation of cells with added RT 01001-CLTH may lead to their apoptosis. Treatment of cells with >5 µM of RT 01001-CLTH in serum depleted media is not recommended, however cells kept in media with serum are resistant to negative effect of RT 01001-CLTH for concentrations above >10 µM.

QUALITY CONTROL: The purity was determined by HPLC analysis. Purity: >95%.

UNPACKING & STORAGE INSTRUCTIONS:

1. Check the container for leakage or breakage.
2. Store at or below -20°C. Repeated freeze/thaw cycles are not recommended.

Note: Please, keep in dark, light may hinder research tool activity.

SAFETY PRECAUTION:

Celther Polska Ltd. highly recommends using protective gloves and clothing. It is also important to avoid breathing vapor, avoid skin contact or swallowing.

WASTE DISPOSAL:

Chemical wastes should always be returned to special company responsible for utilizing such type of wastes.



PRODUCT DATA SHEET

CELther POLSKA SP. Z O.O. WARRANTY

The viability of Celther Polska products is warranted for 30 days from the date of shipment and is valid only if the product is stored according to the information included in this product information sheet. Since this product is under development – we cannot guarantee same functionality in all experimental conditions.

DISCLAIMERS/LEGAL INFORMATION

This product is intended for laboratory research purposes only. It is not intended for use in humans. While Celther Polska uses reasonable efforts to include accurate and up to date information in this product sheet. Celther Polska makes no warranties or representations as to its accuracy. Celther Polska does not warrant that such information has been confirmed to be accurate. This product is sent with the condition that you are responsible for its safe storage, handling, and use. Celther Polska is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of the research tool, Celther Polska is not liable for damages arising from any misidentification or misrepresentation.

While Celther Polska allows to do further modification of the product, Celther does not allow to resale the product without Celther permission. Furthermore, this product is intended to be used in a purchaser laboratory research only.

The above information is believed to be correct, but does not purport to be all inclusive and should be used only as a guide for experienced personnel. Celther Polska shall not be held liable for any damage resulting from the handling or from contact with the above product(s).

CONTACT INFORMATION

All information relating this product is available via email: orders@celther.com

Lodz, Poland, October 2018.